

Reagenzglasgestelle, IDL Test Tube Racks

Reagenzglasgestelle aus rostfreiem 18/10-Stahl, seidenglanz gebürstet, entgratet, punktgeschweißt.

Anzahl Bohrungen	Ø mm	Länge mm	Breite mm	Höhe mm	Best.-Nr.
2 × 10	14	200	48	70	39 00 01223
2 × 12	14	236	46	70	39 00 01234
4 × 12	14	238	84	70	39 00 01244
1 × 5	18	126	32	70	39 00 01201
2 × 6	18	149	55	70	39 00 01212
2 × 10	18	241	55	70	39 00 01225
2 × 12	18	287	55	70	39 00 01237
3 × 12	18	287	78	70	39 00 01241
4 × 10	18	241	101	70	39 00 01231
4 × 12	18	287	101	70	39 00 01246
2 × 6	20	162	60	70	39 00 01216
2 × 10	20	260	60	70	39 00 01227
2 × 12	20	310	60	70	39 00 01239
2 × 6	26	196	72	70	39 00 01218
2 × 6	33	242	86	70	39 00 01219

