

Becher, niedrige Form, DURAN®, DURAN Group / IDL Beakers, Low Form, DURAN®



Becher mit Teilung und Ausguss, DIN 12331, ISO 3819, mit Retrace Code (Chargenzertifikat via Internet).

Volumen ml	A.-Ø mm	Höhe mm	Best.-Nr.	St.	VE St./Pack
5 ¹⁾	22	30	11 21 10607	1-9	10
10 ¹⁾	26	35	11 21 10608	1-9	10
25	34	50	11 21 10614	1-9	10
50	42	60	11 21 10617	1-9	10
100	50	70	11 21 10624	1-9	10
150	60	80	11 21 10629	1-9	10
250	70	95	11 21 10636	1-9	10
400	80	110	11 21 10641	1-9	10
600	90	125	11 21 10648	1-9	10
800	100	135	11 21 10653	1-9	10
1000	105	145	11 21 10654	1-9	10
2000	132	185	11 21 10663	1-9	10
3000	152	210	11 21 10668	1-3	4
5000	170	270	11 21 10673		1
10000 ²⁾	217	350	11 21 10686		1

① Preis pro Stück bei Abnahme kompletter Verpackungseinheiten (VE).

1) Ohne Teilung.

2) Nicht nach DIN und ISO

