

Abdampfschalen, DURAN®, DURAN Group / IDL Evaporating Dishes, DURAN®



WWEQ2044



Abdampfschalen mit Ausguss, DIN 12336.

Rand A.-Ø mm	Höhe mm	Volumen ml	Best.-Nr.	St.	VE St./Pack
50 ¹⁾	25	15	11 21 30132	1-9	10
60 ¹⁾	30	45	11 21 30134	1-9	10
70 ²⁾	35	60	11 21 30138	1-9	10
80 ²⁾	45	90	11 21 30141	1-9	10
95	55	170	11 21 30144	1-9	10
115	65	320	11 21 30149	1-9	10
140	80	600	11 21 30154	1-9	10
190	100	1500	11 21 30159	1-9	10
230	130	2500	11 21 30163	1-9	10

❶ Preis pro Stück bei Abnahme kompletter Verpackungseinheiten (VE).

1) Ohne Aufdruck

2) Ohne Schriftfeld

