

## Korkringe, IDL Cork Rings



WW7R3013



Ringe aus Kork als Unterlagen für Rundkolben, 1a Qualität.

| A.-Ø<br>mm | I.-Ø<br>mm | Höhe<br>mm | für Kolbengröße<br>ml | Best.-Nr.          |
|------------|------------|------------|-----------------------|--------------------|
| 80         | 30         | 30         | 100                   | <b>24 00 20080</b> |
| 110        | 60         | 30         | 250 – 500             | <b>24 00 20110</b> |
| 140        | 90         | 30         | 1000 – 2000           | <b>24 00 20140</b> |
| 170        | 120        | 30         | 2000 – 5000           | <b>24 00 20170</b> |

