

Kühler, Rettberg Condensers

webcode  WWQVD080

rettberg.biz
GLASAPPARATEBAU • LABORATORY GLASSWARE

① Kugelkühler (Allihn-Kühler), Borosilikatglas 3.3, nach DIN 12581

Schliff NS	Mantel-L	Best.-Nr.
mm		
14,5/23 ¹⁾	160	13 40 83120
29/32	250	13 40 83225

Schliff NS	Mantel-L	Best.-Nr.
mm		
29/32 ¹⁾	400	13 40 83240

② Dimrothkühler, Borosilikatglas 3.3, nach DIN 12591

Schliff NS	Mantel-L	Best.-Nr.
mm		
14,5/23 ¹⁾	160	13 40 85116
29/32	250	13 40 85225

Schliff NS	Mantel-L	Best.-Nr.
mm		
29/32	400	13 40 85240

③ Intensivkühler, Borosilikatglas 3.3, nach DIN 12593

Schliff NS	Mantel-L	Best.-Nr.
mm		
14,5/23	160	13 40 86120
29/32	250	13 40 86225

Schliff NS	Mantel-L	Best.-Nr.
mm		
29/32	400	13 40 86240

④ Liebig-Kühler (West-Kühler), Borosilikatglas 3.3, nach DIN 12576

Schliff NS	Mantel-L	Best.-Nr.
mm		
14,5/23	160	13 40 82116
29/32	250	13 40 82225

Schliff NS	Mantel-L	Best.-Nr.
mm		
29/32	400	13 40 82240

¹⁾ nicht nach DIN

