

Laborthermometer, Amarell Laboratory Thermometers

- ① **Einschlussthermometer, geeicht, mit amtlichem Schein**, Kapillare unbelegt, Richter-
verschluss, 8 bis 9,5 mm Ø.

Messbereich °C	Teilung °C	Füllung	Länge mm	Best.-Nr.
-50 – +50	1	Tol	300	16 10 25000
-100 – +30	1	Tol	300	10 10 15010
0 – +100	1	rot	305	10 10 15301
0 – +150	1	rot	305	10 10 15302
0 – +200	1	rot	300	10 10 15303
0 – +50	0,5	rot	220	10 10 15401
0 – +100	0,5	rot	270	10 10 15402



WWERJ004

